

3DISC



User Manual



Software Version 3.8

Document	Software version	Date (d.m.y)
OVO User Manual IOS-FPM-71-001_EN	1.0	01.08.2018
	2.0	01.10.2019
	3.0	01.01.2020
	3.4	20.09.2021
	3.4.x	27.10.2021
	3.4.x	09.11.2021
	3.5	01.02.2022
	3.5.x	08.04.2022
	3.6	15.06.2022
	3.6.x	20.07.2022
	3.7	20.10.2022
	3.7.x	09.11.2022
OvO User Manual	3.8	26.02.2023



The Scan&Tell communication-hub

Introducing OVO	7
1.1 Introduction	7
1.2 3Disc Online Help	7
1.3 Indications supported	7
1.4 Certification and compliance	8
1.5 Regulatory	8
1.6 Intended Use	8
1.7 Classifications	8
1.8 3DISC Privacy Policy and GDPR Privacy Statement	9
2. Components Overview	10
2.1 Hardware Components	10
2.2 Software Components	10
2.3 System Requirements and Specifications	11
2.4 Software output and design software compatibility	12
2.5 Minimum Computer Requirements	12
2.6 Environment Conditions	14
2.7 Power Input	14
2.8 Reusable Tips	15
2.9 Scanner Base and Handpiece	15
2.10 Calibration	15
3. Installing and connecting the OVO	16
3.1 Connecting the OVO	16
3.2 Configuring Windows High Performance Graphics Mode	18
4. Getting Started With 3DiscClinic	20
4.1 Launching 3DiscClinic	20
Registration & Activation	20
Accessing the Admin User Account	21
3 21	
4.2 Setting Up A 3DiscClinic User Account	21
Customizing the Admin User Account	21
4.3 Adding a New User Account	22
Password Protecting a User Account	23
4.4 Linking a User Account to the 3DiscCloud	23
4.5 Accessing the 3DiscClinic Start Screen	26
4.6 Overview of the 3DiscClinic Start Screen	26
Access directly your 3DiscCloud account	27
Enabling Direct Scan	27
5. Configuring Settings & Preferences	31
5.1 Configuring Case Setup Settings	31

5.2 Configuring Lab Sheet Settings	31
5.3 Configuring 3D Settings	31
6. Managing Patients & Cases	32
6.1 Managing Patient Profiles	32
6.2 Managing Patient Cases	32
6.3 Using Direct Scan	34
6.4 Merging a Direct Scan with a Patient Profile	34
6.5 Viewing Case Status	36
6.6 Setting Up a Case in 3DiscClinic	37
7. Scanning with Heron™ IOS	38
Scanning the Maxillary Arch	38
7.1 Using Scan Tools	38
7.2 Using the Quality Map ¶	38
7.3 Taking Live View Screenshots	39
7.4 Using Auto-Realignment ¶	39
7.5 Using 3DiscClinic Live Scan Tools ¶	40
7.6 Carrying out a Manual Bite Alignment	40
8. Finalizing a Patient Case	41
8.1 Using Case Finalization Tools	41
9. Communicating With Labs	42
10. Maintenance	43
9.1 Cleaning the Handpiece	43
9.2 Cleaning and Sterilizing Tips	43
Step-by-Step Procedure	44
9.3 Disposal	45
9.4 Calibration	45
11. Safety Guidelines and Warnings	46
11.1 Warnings and Symbols	46
11.2 General Guidelines	46
11.3 General Warnings	47
System Modification	47
Approved Software	47
Equipment Failure	47
11.4 Mechanical Hazards	47
Moving Parts	47
Dropped Equipment	48
Base	48
11.5 Electrical Safety	48

Electrical Shock	48
Stress on Cables	48
11.6 Eye Safety	49
11.7 Hygiene	49
11.8 Precautions During Systems Operation	49
49	
11.9 Accessories	51
11.10 EMC Guidance and Declaration	51
Electromagnetic Emissions.....	51
Electromagnetic Immunity.....	51

12. Support, Warranty and Repair Service57

12.1 Support	57
Operating Issues Checklist	57
12.2 Standard Warranty	57
12.3 Repair Service	58
12.4 Out of Warranty Repair Service	58



Introducing OVO

1.1 Introduction

Thank you for purchasing the **OVO** from **3DISC**.

The OVO is designed and developed to produce high-quality digital intraoral scans or models, for dental restoration or analysis.

Designed with you as a dentist and your patient in mind, the OVO scanning device is lightweight, versatile and easy to use, enabling fast, accurate scans and enhanced patient experience.

Combined with the 3DiscClinic 3D Scanning and Case Management software, and the 3DiscCloud file-sharing platform for ordering and communicating with labs, the Heron IOS solution provides a comprehensive, intuitive and fully digital experience.

We hope you enjoy your new intraoral scanner solution.

1.2 3Disc Online Help

For your comfort, **3DISC Online Help** is now available at:

3DISC Online Help

Links to 3DISC Online Help are indicated by the following icon:



1.3 Indications supported

The software enables you to select the following indications* when filling out orders for restoration:

- Conventional crowns
- Anatomic crowns
- Copings
- Provisional crowns
- Anatomical pontics
- Reduced pontics
- Provisional pontics
- Inlays/Onlays
- Implant abutments
- Implant-based bridges
- Tooth-based bridges
- Orthodontic aligners
- Nightguards
- Splints
- Retainers
- Bleach trays
- Sleep appliances
- ...

* Verify with your dental lab or service provider about their capabilities to produce different indications.

1.4 Certification and compliance

The system has been tested and conforms to the following standards:

- IEC 60601-1, Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
- IEC 60601-1-2, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances
- Requirements and tests:
- IEC 62471, Photobiological safety of lamps and lamp systems

1.5 Regulatory

The OVO is manufactured and marketed in accordance with US FDA Regulations and EU Medical Device Regulation 2017/745.

1.6 Intended Use

The OVO is an optical impression system. It is used to record the topographical characteristics of teeth, dental impressions, or stone models for use in the computer-aided design (CAD) and computer-aided manufacturing (CAM) of dental restorative prosthetic devices.



WARNING: Unintended use of the system can result in physical injury to the patient and operator, and damage to the system.



CAUTION: Rx only - Federal law restricts this device to be sold by or on the order of a Dentist.

1.7 Classifications

The OVO system has the following classifications

- Protection against electrical shock: Type B Applied Part
- Protection against harmful ingress of water: Ordinary equipment (IPX0)
- Safety of application in the presence of a flammable anesthetic material with air or with oxygen or nitrous oxide: Equipment not suitable for use in

the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

1.8 3DISC Privacy Policy and GDPR Privacy Statement

Your privacy and the privacy of the patient data managed by you is important to 3DISC.

For information on 3DISC Privacy Policy, go to:

[Privacy policy - 3DISC](#)

To access 3DISC's GDPR Privacy Statement, go to:

[GDPR Privacy statement - 3DISC](#)


2. Components Overview

The OVO system is composed of hardware and software components.

2.1 Hardware Components

The Heron IOS system is composed of the following hardware components*:

1. Heron IOS Scanner
2. Heron IOS base
3. Removable autoclavable tips
4. USB 3.0 cable
5. AC/DC power adapter



Hardware Component	Part Number
OVO (Includes 3 Tips)	IOS-FP-71-001
USB 3.0 Cable	IOS-CP-00-043
AC/DC Power Adapter	IOS-CP-00-088

Save the Box: It is highly recommended that you store the packaging in a safe place and do not dispose of it. The original packaging box is optimal for any necessary transportation or shipment of the OVO.

***Note:** *Packaging details and contents may vary from those described in this guide.*

2.2 Software Components

The OVO system comprises the following software components:

- 3DiscClinic™ : 3D Scanning and Case Management software.
- 3DiscCloud™ : Dedicated cloud platform for ordering and communicating with labs.

2.3 System Requirements and Specifications

Scanner Specifications

Scanner Type	Hand-held (chairside) scanner that creates optical impressions for dental restorations.
Design	Compact, lightweight, ergonomic – designed to be operated with little physical effort.
Dimensions Base	Size: L 306mm, W 98mm, H 72mm
Dimensions Scanner	Weight: 150 grams
	Size: L 256mm, W 43mm, H 43mm
	Cable length (scanner to base): 2m
Power Requirement	DC 5.0V / 4A (Power supply included)
Scanner Tip	Reusable up to 250 times, sterilize using steam autoclave
Heating Element	Ventilated. Prevents formation of fog on optics
Acquisition Method/Imaging technology	Hybrid technology: active stereo imaging and structured light
Sensor technology	CMOS
Color Scanning	24-bit (8-bit per channel)
Scanning frequency	25-30 FPS
Imaging field-of-view	12mm x 14mm
Light sources	High-power LEDs

Scanning Process

Tooth Preparation	No powder or spray required
Scanning Principle	Continuously scanning and accumulating (stitching) depth and color data
Distance Scanner - Tooth	-1mm – 19mm
Possible contact duration by operator	<10 min. Note: May vary with hardware configuration
Operator accessible part	Handpiece
Possible contact duration by patient	$t \leq 10$ min

Patient accessible part (Type B Applied Part)	Tip (autoclavable)
Computer – Scanner Interface	USB 3.0

2.4 Software output and design software compatibility

Output File Format	STL, PLY, OBJ
Compatibility with CAD/CAM Systems	Open Architecture Output format STL, PLY, OBJ Compatible with most Dental CAD systems

2.5 Minimum Computer Requirements

The following requirements have been defined by 3DISC to ensure the 3DiscClinic software operates properly. The related configurations have been tested by 3DISC.

For an online version, visit: [Minimum Computer Requirements](#)



IMPORTANT NOTICE: The use of any other hardware units and/or other base software modules to run the 3DiscClinic software is not recommended and is not supported by 3DISC.

Minimum Software Requirements

Operating System	Windows 10 (Excluding Windows 10 S, now defunct) Administrative rights required.
Disk Space	100 GB or greater of free disk space
Ports	At least 1 x USB 3.0 port (SuperSpeed)
Nvidia Driver	Nvidia Studio driver version 516.94 is currently required. Nvidia gamer-ready driver should not be used with the HeronClinic software.
NVIDIA GPU operation mode	The PC must be able to ensure that the Nvidia GPU is the only GPU activated on the PC and that <u>any Intel integrated GPU is disabled</u> . (This is usually achieved using high level configuration tools or BIOS commands). Warning: Certain models of computer of the following brands - Acer, MSI, HP,... - provide no option to only have the Nvidia GPU active. This may affect performance: low FPS or freeze during scan.

Screen resolution	Full HD (1920 x 1080) with DPI 100% NOTE: The use of 4K (3840 x 2160 pixels) or Ultra-Wide (3440 x 1440 pixels) resolutions is also possible, however the impact on performance has not been quantified by 3DISC at this stage.
-------------------	--

Software Configuration Recommendations

Windows automatic updates	3DISC recommends deactivating all Windows automatic updates (except for security updates).
Nvidia driver automatic updates	Nvidia driver automatic updates <u>should be disabled</u> .
Windows Battery Settings	On laptops, the battery setting in Windows should be configured to high performance mode only , with <u>no battery saving option</u> .

Minimum Hardware Requirements

CPU Type	Intel 10, 11 and 12 generation. Intel i7 or i9 – 4 Cores give best performance.
CPU Clock	2.8 GHz clock or greater “Turbo” and “boost” speeds cannot be considered.
Memory	32 GB of RAM or greater (DDR4 or better)
Graphics Card Memory	6GB of RAM minimum on the graphics card are needed. Below this, the software will not launch: an error message will inform you that the minimum requirement is not reached.
GPU	<ul style="list-style-type: none"> • Quadro RTX3000, RTX4000 and above for laptop and desktop • RTX2070 for laptop and desktop • RTX2080 for laptop and desktop • RTX2080TI for desktop • RTX3070 for laptop and desktop • RTX3080 for laptop and desktop • RTX3090 for desktop

These PC requirements may be revised without notice by 3DISC to take into account observations made on the field or additional test results performed by our teams.



IMPORTANT: Compatibility of AMD GPUs is not guaranteed with the Heron™ IOS.

Not meeting minimum hardware requirements will affect the performance of the scanner.

2.6 Environment Conditions

Operating Temperature	10°C to 30°C
Operating Relative Humidity	10% to 80% (non-condensing)
Storage Temperature	- 20°C to 60°C
Storage Relative Humidity	10% to 80% (non-condensing) Indoor use only
Installation Category	1
Pollution Degree	2
Ingress of Liquids	IPX0
Protective Class	Class IIIb
Overvoltage category	II per IEC 60664-1
Max. working condition	Continuous cycles with image capture and transmission from/to Notebook or non-medical grade PC.
Other possible accessories (IEC60601-1 3rd, Cl. 16)	Notebook with AC/DC Adapter.
Equipment Maintenance	No user maintenance is required, and no user service is allowed. Please contact technical support in case of problem.
Cleaning	Do not try to clean the inside of the device. Refer to section 9 : Cleaning the Handpiece for cleaning and sterilization.

2.7 Power Input

The power adapter input is 5V DC, 100-240V AC, 50-60Hz.

2.8 Reusable Tips

Scanner tip is autoclavable up to 250 times in a steam autoclave when used with min cycle:

- 132°C (270°F) at 4 minutes, or
- 134°C (273°F) at 4 minutes, or
- 121°C (250°F) at 45 minutes.

See below: [section 9.2 Cleaning and Sterilizing Tips](#).

2.9 Scanner Base and Handpiece


The scanner body consists of the Docking Base and Handpiece, which are connected by a flexible, non-detachable cable.

2.10 Calibration

The OVO™ is factory calibrated. In the case of calibration issues due to transport, please contact your reseller or 3Disc support technician.

3. Installing and connecting the OVO

3.1 Connecting the OVO



For an online version, visit: [Getting Started with OVO](#)



To install and connect the OVO Scanner:

- Step 1.** Place the docking base on a flat, stable surface and place the OVO handpiece securely on the base.
- Step 2.** Connect the AC/DC power adapter cable to the docking base (the connector socket is located underneath the base of the scanner).



WARNING: Make sure you use the 5.0V 4A power adapter provided. Failure to do so may result in damage to the scanning device.

- Step 3.** Connect the provided USB 3.0 cable to the docking base (the connector socket is located



WARNING: Using a USB cable other than the one provided may result in system malfunction or reduced performance.

underneath the base of the scanner).

Step 4. Connect the other end of the USB 3.0 cable to the computer.



IMPORTANT: Make sure to use a USB port that is compatible with USB 3.0 (SuperSpeed), usually indicated by this symbol: not doing so may result in system malfunction or reduced performance.



IMPORTANT: When using a desktop computer, it is strongly recommended to plug the USB cable to a USB port located at the back of the computer; not doing so may result in system malfunction or reduced performance.

DO NOT plug the USB cable into an intermediate hub.

Step 5. Connect the adapter block provided to a power outlet.

IMPORTANT: Connect to Power Supply before Scanning!

If your laptop computer is not properly connected to a power outlet, the following message will appear:

In this case, please connect your laptop to the power supply **before proceeding to use to scanner.**



This system is not properly connected to AC Power or is running on battery. This will cause the software to run poorly or very slowly.

Please connect the power supply to this system before proceeding to use this software.


OK



IMPORTANT: If using a laptop computer, make sure the power supply is connected to a power outlet and not running on battery power. Failure to do so will mean that the scanner will not have sufficient power to produce images.

3.2 Configuring Windows High Performance Graphics Mode

On laptop computers, battery settings in Windows should be configured to high performance mode only, with no battery saving option (Settings/System/Battery).



To select High Performance mode for the HeronClinic application:

1. In Windows **Graphic settings**: click to activate **Hardware-accelerated GPU scheduling**
2. In **Graphics performance preference**: select **Desktop app**
3. Click **Browse** and select the **3DiscClinic** app as shown (**Launcher.exe**):



4. In **Options**, select **High performance**, and click **Save**.
5. Repeat this for the **HeronClinic PatientDB.exe** and **Scan3D.exe** executable files.
6. Restart your PC to apply changes.



IMPORTANT: On laptop computers, battery settings in Windows should be configured to high performance mode only, with no battery saving option (Settings/System/Battery).

WARNING: Unsuitable installation sites



WARNING: Unsuitable installation sites:

- Locations with excessive humidity or dust
- Locations subject to high temperature
- Locations subject to shaking or vibration
- Locations exposed to considerable electrical or magnetic noise, or other forms of electromagnetic energy

4. Getting Started With 3DiscClinic

For an online version, visit: [Getting Started with OVO](#)



4.1 Launching 3DiscClinic

Once you have connected the OVO scanner, you are ready to launch the 3DiscClinic software installed on your computer.

Step 1. Click on the **3DiscClinic** desktop icon to launch the 3DiscClinic software.



Registration & Activation

On first launch, you will be invited to register and activate the 3DiscClinic software on your computer.



Device Registration

* First Name

* Last Name

* Email Address

Scanner S/N 10426

Installation Date 5/23/2022

* Country United States

* Address

Address (continued)

* City

State/Zip Code

* Phone +1

By using this form, you agree to the storage and processing of your data.


SUBMIT

In the Device **Registration** dialog box:

Step 2. Enter your registration details (required information is indicated by *)

Step 3. Tick the **consent** checkbox, to consent to the collection and processing of data.

Step 4. Click **SUBMIT**.



IMPORTANT: It will not be possible to perform new scans or export existing scans if the device has not been registered on the PC.

Data collected is used in accordance with the **General Data Protection Regulation (GDPR)** and is not shared with third parties. For information, see:

- [section 1.8 - 3DISC Privacy Policy and GDPR Privacy Statement](#)
- <https://gdpr-info.eu/>

Accessing the Admin User Account

The 3DiscClinic Start Screen will display the default 3DiscClinic Admin User Account (**HC**) in the left-hand menu.

To access the Admin user account:



- Step 1.** Click on the Admin User Account (**HC**) icon in the left-hand menu.



The HeronClinic Sign In page invites you to **sign in** or **create a new account**.

- Step 2.** Enter your password and click **Sign In**.

The Admin User Start Screen displays as shown:



4.2 Setting Up A 3DiscClinic User Account


Customizing the Admin User Account

To customize the **3DiscClinic** Admin User Account (**HC**):

- Step 3.** Click the System icon located in the top right-hand corner of the User Start Screen, to access the **Settings** interface.





Step 4. Select **Users** in the left-hand **Settings** menu, and click on **Add and Edit Users**.



Step 5. Select the default “3Disc Clinic” User profile, and click on the **Edit** icon to customize the default Admin User account.

Note: The default admin ID “3DiscClinic” cannot be modified.

Step 6. In the **User Image** field, you can click to open the **Picture** dialog box, to add or take a User Photo.



Step 7. To apply changes, restart the **3DiscClinic** application.

4.3 Adding a New User Account

Only the **3DiscClinic** Admin User profile (HC) can add new User Accounts.




To add a new User account:

Step 1. In Settings/Users, click **ADD NEW USER**:

This opens the **New User** dialog box.



Step 2. Enter User information: email, first name, last name, registration no. (optional), photo (optional).



Password Protecting a User Account

To password protect a User profile (recommended):

- Step 3.** Select **Yes** in the **Use Password** field.
- Step 4.** **Enter and confirm** the password to apply to this user account.
- Step 5.** Click **CREATE**.
- Step 6.** To apply changes, **close and restart** the application.



IMPORTANT: To protect the privacy of patient data processed by you, it is important to password protect all user accounts.

4.4 Linking a User Account to the 3DiscCloud

You can link your **3DiscClinic** User account(s) to the **OVO** solution's dedicated **3DiscCloud** platform. The **3DiscCloud** enables users to easily and efficiently setup and manage file transfers and connections with labs.



Note: The Heron IOS **3DiscCloud** platform is designed to facilitate file sharing and manage connections with laboratories. **It is not a cloud storage service.**

To link the **3DiscClinic** User profile to a **3DiscCloud** account:

- Step 1.** Using the Admin account, In Add and Edit Users, click **ADD NEW USER** to create a User account.
- Step 2.** Enter account details, and in the **Use HeronCloud** field, Select **Yes**.
- Step 3.** Click Create, to add the new User account.
- Step 4.** Click **OK**.

The screenshot shows a web form titled "New Account" with the HeronCloud logo at the top. The form contains the following fields: Email (with a dropdown arrow), Password, Repeat Password, First Name, User, Last Name, Name, Role (with a dropdown arrow), Institution, Country (with a dropdown arrow), Phone, and Address. A red "Create" button is located at the bottom right of the form.

The **3DiscCloud Login** window will open automatically, inviting you to **Login** or to **Create a New 3DiscCloud Account**.


For more information, see below:

Chapter 0 -

Communicating With Labs:

- **Error! Reference source not found.** - Error! Reference source not found.
- **Error! Reference source not found.** - Error! Reference source not found.

4.5 Accessing the 3DiscClinic Start Screen




User accounts are displayed in the left-hand menu of the 3DiscClinic Home Page.



To access the User **Start Screen**:

Step 1. Click on your User profile in the left-hand menu of the 3DiscClinic Home page.


The **3DiscClinic Login** window will open automatically, inviting you to **Login** or to **Create a New Account**.



Step 2. Enter your password and click **SIGN IN**.

4.6 Overview of the 3DiscClinic Start Screen

When you login to your **3DiscClinic** User account, the Start Screen is displayed as shown:



A. Patients/Cases View

You can toggle between **Patients** and **Cases** views, and Filter or search desired Patients/Cases.


Filter/Search

Search by term or Click on the Filter icon to open the **Filters** dialog box.



You can display a list of Patients by:

- **Gender**
- **Case Creation Date** (or range of dates)
- **Date of Birth**



Select A User Start Page Default View

To select a preferred **Default View** for the **User Start Page**:



Step 1. Click **Settings (E)**.


Step 2. In the left-hand **Settings** menu select **General**.

Step 3. Click **Preferred Start Page**, and

Step 4. In the drop-down menu, select **Recent Patients** or **Recent Cases**.

Access directly your 3DiscCloud account

You can access the 3DiscCloud dashboard linked to your account by clicking on the Cloud icon opposite your username.



Enabling Direct Scan

It is possible to directly access the scanning module without first selecting a patient profile or adding a case.


To enable **Direct Scan**:




Step 1. Click **Settings (E)** to access System Settings.

Step 1. In the Settings Menu, select **General**.

Step 2. In General settings, Select **Enable Direct Scan Workflow**. (Note: this option is disabled by default).



The **Start Scanning** bar is displayed in the User Start Screen.



Next Steps

For information on configuring **3DiscClinic Settings**, see:

- [Chapter 5 - Configuring Settings & Preferences](#)

For information on managing patient cases in **3DiscClinic** see:

- [Chapter 0 -](#)

- [Managing Patients & Cases](#)

For information on scanning with **3DiscClinic** see:

- [Chapter 7 -](#)

- Scanning with HeronTM IOS

5. Configuring Settings & Preferences

For an online version, visit: [HeronClinic Settings](#)



5.1 Configuring Case Setup Settings

For an online version, visit: [Case Setup Settings](#)



5.2 Configuring Lab Sheet Settings

Lab Sheet settings enable users to customize Order forms for use by their organization, by adding the name, address and corporate logo to their Lab Sheets, and selecting a page format (A4/US Letter).

The screenshot shows a form titled "Lab Sheet" with the following fields:

- Institution Name
- Institution Address
- Page Size: A4 (with a dropdown arrow)
- Institution Logo (with a dropdown arrow)

Order Forms are published in PDF format, and created automatically based on the **Case Setup Options** and **Lab Sheet** settings selected by you.

5.3 Configuring 3D Settings

For an online version, visit: [Configuring 3D Settings](#)



6. Managing Patients & Cases

For an online version, visit: [Managing Patients & Cases](#)



6.1 Managing Patient Profiles

For an online version, visit: [Managing Patients](#)



6.2 Managing Patient Cases


Previewing a Patient Case

The **Case Finalization page** provides detailed case information at-a-glance: creation date, case ID, restoration types, status, and any 3D models and 2D images associated with the case...

To preview a patient case:

Step 1. Select a patient in the left-hand menu.

Step 2. Select a case in the list of cases associated with the patient profile (4).



The **Case Preview** is displayed as shown.

Opening a Patient Case

To open an existing patient case:


Step 3. Click **Open Case** in the **Case Preview** window.



This will open the **Finalization** page.

To directly access any step of the **Case Workflow**:

Step 4. Click on the relevant icon in the left-hand




menu.

For information on using the **Finalization** page, see below:

- [Chapter 41 - Finalizing a Patient Case](#)

6.3 Using Direct Scan

It is possible to skip the **Case Setup** step and directly access the Heron™ IOS Scanning module. To do so, select the **Start Scanning** click bar in the 3DiscClinic Start Screen.



In the case of a Direct Scan, the scan data is not associated with an existing patient profile, and is attributed an automatically generated patient ID.

Note: *The Direct Scan feature is not activated by default. To activate click **Enable Direct Scan Workflow** in System Settings.*

Enable Direct Scan Workflow



6.4 Merging a Direct Scan with a Patient Profile

It may be useful to merge a direct scan you have carried out with an existing patient profile. You do this by merging the automatically generated DirectScan profile with an existing profile.

To merge patient profiles:


Step 1. Select the profile to be merged and click the icon to the right of the patient folder (see above).

Step 2. Select merge in the expanded list.

In the **Merge Patient** window:

Step 3. Search and select a patient ID.

Step 4. Click **Merge** and **Confirm**.




Cloning a Case

It is possible to duplicate an existing case and all associated data, including order form, scan data, intraoral photos, and attachments.

To clone a case:

Step 1. Click on the 3-point icon to the right of the selected case as shown



Step 2. Click **Clone** in the drop-down list.

Note: You cannot clone a case for which scan data has already been deleted.

Moving a Case

If a case has been associated with an incorrect patient profile, you can move a case from one patient profile to another.

To move a case:

Step 1. Click on the 3-point icon to the right of the selected case.


Step 2. Click **Move** in the expanded list.

In the **Move Case** dialog box:

Step 3. Select or create a patient profile (**ADD NEW PATIENT**).

Step 4. click **MOVE**.

The case is associated with the selected patient profile.



6.5 Viewing Case Status





To access **Case status** details in the **Finalization** page:



Step 1. Click on the status icon as shown above.

The expanded list provides information on the User, Case ID, Lab, Date the case was last modified, and the send status.

Case Status Icons

Case Status	Description
 Created	Case is created and currently in progress.
 Scanned	Case is scanned and finalized.
 Uploaded	Case is uploaded to lab via 3DiscCloud.
 Closed	Case is closed. When the Auto Clean Raw Scan Data option is enabled in Data Settings , cases for which raw scan data is deleted are automatically closed.

6.6 Setting Up a Case in 3DiscClinic

The following workflow briefly describes the typical steps and options available to set up a patient case in **3DiscClinic**.

For an online version, visit: [Creating a Patient Case](#)



For an online version, visit: [Case Setup Settings](#)



7. Scanning with Heron™ IOS

For an online version, visit: [Scanning with Heron IOS](#)



Scanning the Maxillary Arch

For an online version, visit: [Scanning the Maxillary Arch](#)

When scanning the Maxillary arch, we recommend you scan in the following order:



7.1 Using Scan Tools

At each step of the scan workflow, you can use the features available in the right-hand **Scan Tools** menu.

For an online version, visit: [Using Scan Tools](#)



7.2 Using the Quality Map ¶

The **Quality Map** tool enables you to assess in real-time if enough data has been collected in the area of interest of the scan.

For an online version, visit: [Using the Quality Map®](#)




7.3 Taking Live View Screenshots

To obtain Live View Screenshots during the scan:

1. Press the **C** key on your keyboard or click on the icon in the scan tools menu.



2D Live View Screenshots are displayed in thumbnail format in the bottom left-hand corner of the screen.



2. Click on a *thumbnail image* to inspect the 2D images in **fullscreen** format.
3. Click on the *Delete* icon to delete as necessary.
4. Click on the *Close* icon in the top right-hand corner to close the fullscreen view.

Live View Screenshots taken during the scan can be viewed in the **Finalization** page.

7.4 Using Auto-Realignment ¶

The **Auto-Realignment** tool carries out realignment of the scans, optimizing the scan data based on **Artificial Intelligence (AI)**. To simplify, the scan data is cleaned by discarding frames that are calculated to be misaligned with respect to the totality of the available data.



You are advised to use the **Auto-Realignment** tool before proceeding to the **Bite Alignment** step, as optimization may modify the occlusal position of the realigned arches.

7.5 Using 3DiscClinic Live Scan Tools ¶

During the **3DiscClinic Scan Workflow** , **Live Scan Tools** are displayed below the digital 3D Model.

For an online version, visit: [Using 3DiscClinic Live Scan Tools](#)

For an online version, click: [Automatic Bite Alignment](#)



7.6 Carrying out a Manual Bite Alignment

For an online version, click: [Manual Bite Alignment](#)



8. Finalizing a Patient Case



For an online version, click: [Finalizing a Case in HeronClinic](#).



8.1 Using Case Finalization Tools

For an online version, click: [Using Case Finalization Tools](#)



9. Communicating With Labs

For an online version, visit: [Communicating with Labs](#)



10. Maintenance

9.1 Cleaning the Handpiece

The entire body, cord and base of the scanner must be wiped down using a **Federal Environmental Protection Agency (EPA)** approved disinfectant that is labeled and specified for tuberculocidal/ mycobactericidal activity. Do not use disinfectant on the nozzle.



IMPORTANT: All components of the scanner (excluding the tips) must be **wiped down** and not sprayed. Avoid getting any moisture, alcohol or disinfectant inside the open scanner chamber.

Recommended and approved surface disinfectants:

- Birex Wipes: TB Claim= 10 minutes- ‘Phenolic (Dual) Water-Based’
- Prospray Wipes: TB Claim=10 minutes- ‘Phenolic (Dual) Water-Based’
- Cavicide Wipes: TB Claim=3 minutes- ‘Phenolics (Dual) Alcohol-Based’
- DisCide ULTRA Wipes: TB Claim-1 minute- ‘Phenolics (Dual) Alcohol-Based’
- Maxiwipe Germicidal Cloth: TB Claim=5 minutes- ‘Phenolics (Dual) Alcohol-Based’
- Ster 1 Plus: TB Claim=5 minutes- ‘Quaternary ammonium and Alcohol-Based’

9.2 Cleaning and Sterilizing Tips



IMPORTANT: The included tips must be autoclaved prior to use as they do not come pre-sterilized.



IMPORTANT: Ensure that the surface of the mirror does not show residues, smudges, scratches, or any damage, as this would affect the performance of the device.

Step-by-Step Procedure


Step 1. Clean the tip with soapy water, ensuring that the mirror is clean and free of smudges, stains, or any residue. Avoid using abrasive cloth materials as this will scratch the mirror.



Step 2. After drying the tip exterior, carefully dry the interior and mirror with lint free wipes ensuring you do not scratch the surface. The mirror should be free of any noticeable debris or water spots.



Step 3. Insert and seal the tip into a sterilization pouch. Make sure the seal is airtight. Each tip should be packaged individually.



Step 4. Sterilize the wrapped tip in a steam autoclave at the following parameters:

- 132°C (270°F) at 4 minutes, or
- 134°C (273°F) at 4 minutes, or
- 121°C (250°F) at 45 minutes



Step 5. Ensure the dry cycle is complete prior to removing the tip from the autoclave. If the pouch is damp with moisture, proper sterilization cannot be guaranteed.



WARNING: Always autoclave the tip wrapped in a sealed sterilization pouch; failure to do so will result in permanent stains on the mirror.



NOTE: Tips should not be placed in an ultrasonic cleaner or any cold sterile solutions. The sterilant solutions will leave a sticky residue or film on the mirror when drying.



WARNING: Do not autoclave the handpiece of the device.



WARNING: Do not remove the pouch before the sterilizer completes its full dry cycle. If the pouch is wet or has any signs of moisture, this can potentially leave water spots on the mirror which can affect image quality during scanning.

9.3 Disposal

The Heron is an electrical device with electronic components inside and should be disposed of in accordance with local environmental laws and regulations.

9.4 Calibration

The OVO is calibrated in the factory and therefore does not require calibration when installed.



WARNING: General prohibition indication. The functionality of the system can be destroyed in the case of incorrect use. If unauthorized changes have been made to the delivered system and accessories, the warranty by 3DISC becomes void. 3DISC will not accept any responsibility or liability for the improper functioning of the product in such a case.



WARNING: Use extreme caution when cleaning the mirror as it is very delicate and is prone to scratching.

If the **OVO** begins to have problems scanning and recognizing teeth models, contact your dealer or **3DISC Support** technician.

If the scanner cannot be recalibrated remotely, this may result in the system being returned for repair/ calibration.

For more information, see below:

- [Chapter 12 - Support, Warranty and Repair Service.](#)

11. Safety Guidelines and Warnings

11.1 Warnings and Symbols



NOTE: Notes represent information that is important to know but which do not affect the functionality of the system.



WARNING: The functionality of the system will be limited in the case of incorrect use.

11.2 General Guidelines

- Do not spill liquids on the body of the device
- Never operate the device in a wet environment.
- Keep the device away from radiators and heat sources.
- Use the device only with the accessories supplied.
- Do not alter the device or open enclosures.



WARNING: General prohibition indication. The functionality of the system can be destroyed in the case of incorrect use. If unauthorized changes have been made to the delivered system and accessories, the warranty by 3DISC becomes void. 3DISC will not accept any responsibility or liability for the improper functioning of the product in such a case.

If any of the following conditions occur, unplug the device from the electrical outlet and contact authorized service personnel:

- The power cord or power adapter is damaged.
- The device has been exposed to water.
- The device has been damaged.
- The device does not operate correctly when the operating instructions are followed.

11.3 General Warnings

System Modification



WARNING: Modifying the system may result in physical injury to the patient and operator, and damage to the system.

Approved Software

The OVO device is designed to operate with the 3DiscClinic software.



WARNING: The Heron™ IOS scanner should only be used with approved, compatible software.

Equipment Failure



WARNING: In case of **system** malfunction or failure, you should:

- Prevent any contact between the system and the patient.
- Unplug the system from the power outlet and the computer.
- Store the system away so it cannot be used by someone else.
- Contact service personnel.

11.4 Mechanical Hazards

Moving Parts



NOTE: All moving parts are inside handheld scanner so do not open the unit.

Dropped Equipment



WARNING: If the scanner tip is dropped, ensure that the mirror is not damaged and that it is not detached; if the tip is damaged it should be disposed of immediately. If the scanner handpiece is dropped or bumped, ensure that no part of the system is damaged as it could affect performance.

Base



NOTE: When not in use, always rest the handpiece on the Base. The Base may be mounted on the wall per provided instructions. Do not place the Base on a slanted surface. Place the cables (power cable and USB cable) where people cannot accidentally get caught in them and potentially damage the system.

11.5 Electrical Safety

Electrical Shock



WARNING: There is a risk of electrical shock when opening or attempting to open any part of the system; only qualified service personal should open parts of the system.

Stress on Cables



WARNING: Care should be taken not to apply unnecessary stress on the cables of the system, whether it is the power cable, the USB cable or the cable between the handpiece and the Base.



WARNING: Only use the power adapter supplied as a part of the system.

11.6 Eye Safety



WARNING: During operation, the system emits a bright, flashing light from its tip. Although the system complies with standard IEC 62471 (Photobiological safety of lamps and lamp systems), prolonged exposure to flashing light may result in discomfort, seizure or eye irritation.

11.7 Hygiene



WARNING: In order to maintain safety for the patient, wear surgical gloves when handling any parts of the system. Always ensure that the tip is mounted on the handpiece before inserting it into the mouth of the patient. Before using the system with a new patient, ensure that the system is disinfected, and the tip is sterilized.

11.8 Precautions During Systems Operation



WARNING: The Heron™ IOS system contains delicate optical and mechanical elements and therefore should be handled carefully. Do not drop, bump, or shake the handpiece or the tip. Always place the handpiece on the Base when not in use. Do not put stress on the cable connecting the handpiece to the Base. Do not submerge the handpiece or the Base in any liquid. Do not place the handpiece or the Base on wet or heated surfaces. Hold the handpiece with a firm grip when handling it.



WARNING: Portable RF communication equipment (including peripherals such as antenna cables and external antennas) should not be used within 30 cm of any part of the Heron™ IOS, including cables specified by the MANUFACTURER. Otherwise, performance degradation of this equipment may occur.



NOTE: During operation of the system, the handpiece and the tip may get slightly warm; this is normal.



WARNING: In order to prevent over-heating of the system, the ventilation opening at the bottom of the handpiece should never be obstructed.



WARNING: Use of this equipment adjacent to or on other equipment must be avoided as it may result in improper operation. If this use is necessary, it is advisable that this and the other equipment be observed to verify that they are operating normally.

OVO Performance

Under normal use, the OVO should transmit images to the laptop/workstation with the installed image manipulation software and the video stream is visible in the bottom left corner of the screen. If not used properly, there may be loss of transmission of image information or slow transmission temporarily.

11.9 Accessories



WARNING: Only use 3DISC approved accessories. Not using 3DISC approved accessories may result in deterioration of performance.

Accessory	3DISC Part Number
USB 3.0 Cable	IOS-CP-00-043
AC/DC Power Adapter	IOS-CP-00-088

OVO is intended for use in a professional healthcare setting with electromagnetic environment specified below.

11.10 EMC Guidance and Declaration


Electromagnetic Emissions

Emissions Test	Compliance	Electromagnetic Environment Guidance
RF Emissions CISPR 11	Group 1, Class A	OVO uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Complies	OVO is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Electromagnetic Immunity

Immunity Test	Test Level	Compliance Level	Electromagnetic Environment Guidance
---------------	------------	------------------	--------------------------------------

Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV /Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV /air	± 8 kV /Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV /air	Surface should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical professional healthcare facility environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	Mains power quality should be that of a typical professional healthcare facility environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0,5 cycle 40% UT (60% dip in UT) for 5 cycles 70 % UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	<5% UT (>95% dip in UT) for 0,5 cycle 40% UT (60% dip in UT) for 5 cycles 70 % UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	Mains power quality should be that of a typical professional healthcare facility environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a location in a typical professional healthcare facility environment.
Radiated RF IEC 61000-4-3	Table 9 in IEC-60601-1-2 2014	Table 9 in IEC-60601-1-2 2014	Portable and mobile RF communications equipment should be used no closer to any part of the OVO system, including cables, than the recommended
Conducted RF IEC 61000-4-6	3V 0.15-80MHz	3V 0.15-80MHz	

	<p>6V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1KHz</p>	<p>6V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1KHz</p>	<p>separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:</p> $d = \left[\frac{3.5}{V^2} \right] \sqrt{P} \quad 150 \text{ kHz to } 80 \text{ MHz}$ $d = \left[\frac{3.5}{E^2} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E^2} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--	---	---	--

NOTE: UT is the A.C. mains voltage prior to application of the test level.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the OVO System that is not Life-Supporting

OVO System is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the OVO system as recommended below, according to the maximum output power of the communications equipment.




















Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{V1}\right]\sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{E1}\right]\sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E1}\right]\sqrt{P}$
0,01	0.12	0.12	0.23
0,01	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

11.11 Symbols on the Heron IOS and Transport Case

Symbol	Description
	Manufacturer's trade name and address (ISO 15223-1)
	Date of manufacture (ISO 15223-1)
	Equipment Power On/Off (push/push)
	USB 3.0 plug
	Warning, Consult Accompanying Documents
	General mandatory action manual
	General prohibition indication
	User Manual Reference
	Directive on Waste Electrical and Electronic Equipment
	Authorized Representative in the European Community
	Warning label for LED
	Non-ionizing electromagnetic radiation
	Direct Current
	Type(B) Level of protection against electric shock
	Consult operating instruction for use.
	European Conformity mark
	Prescription symbol
	Non-sterile (Scanner-Tips) (IOS-FPL-71-001)
	User manuals are available electronically at the link provided (3disc.com/support-resources/heron-ios-user-manuals/)



INMETRO Certification Mark



Ukraine Conformity Mark

12. Support, Warranty and Repair Service

12.1 Support

If you have questions about the software, please consult the manual and Help menu in the software. If you are experiencing issues with your software, please check the list of common issues provided below prior to contacting a dealer. It could be simply a question of a minor issue that can be fixed quickly. However, if you're still experiencing problems after following the recommendations in this section, then please contact the dealer where you bought the equipment.

Operating Issues Checklist

Issue	Recommendation
There's a memory full error message that pops up when the software is open.	Clear some space on the C Drive
The status in the Live view window is "Disconnected".	Check that you have external power to the Heron and that the USB cable is connected to a USB 3 Port.
Scanning is very slow.	Check that the Laptop is connected to an external power source.
The corners are cut in the live view window.	Check that the Tip is correctly mounted and when rotating it is locking in place with a click.
There is a red square in the scan window	Go back to a tooth that is scanned and start from there again
No images appear when scanned but everything else (e.g. live window image, sounds, FPS) works fine.	The scanner might need to be recalibrated. Please contact your local dealer for support.
There are spots on the Live view window.	Check and clean the mirror of the tip.
Where can I get the OVO software and manuals?	Software and manuals can be downloaded in the Support section of the 3DISC website.

12.2 Standard Warranty

3DISC warrants its non-consumable hardware products to be free from defects in materials and workmanship. The warranty covers the cost of parts and labor to repair the product.

Please keep the shipping container for future use. Products returned to the factory for

repair should be properly packaged. To obtain warranty service, follow the procedure described in the Repair Service section. Failure to do so will cause delays and additional expense to the customer.

The warranty is valid when the product is used for its intended purpose and does not cover products which have been modified without written permission from 3D Imaging and Simulation Corp. Americas, or which have been damaged by abuse, accident or connection to incompatible equipment.

This warranty is in lieu of all other warranties, expressed or implied.

12.3 Repair Service

The OVO cannot be serviced locally. In the event of a hardware malfunction, contact your dealer to arrange for a swap unit (same model or newer) so your unit can be replaced, and work can continue. Some testing might be needed in order to verify the Hardware/Software error or malfunction.

The company reserves the right to cease providing repair, maintenance, parts and technical support for its non-consumable hardware products five years after a product is discontinued.

12.4 Out of Warranty Repair Service

Out of warranty repair service is available in selected geographical locations. Contact the supplier for current terms and rates.

We hope this User Manual was helpful to you.
For additional material and user information go to:

3disc.com/support-resources/

OVO User Manuals
OVO How-to Videos
OVO Training Videos

3DISC



**3D Imaging and Simulations Corp.
Americas**

365 Herndon Pkwy #18
Herndon, VA 20170
USA

Tel: +1 703 430 6080

Tel: +1 800 670 0363 (Toll free)



3DISC Europe

3DISC Dental Connect

191, avenue Charles de Gaulle,
92200, Neuilly-sur-Seine,
France

Tel : +33 (0)1 42 25 73 98

Contact 3DISC

Info@3DISC.com

Contact support

support@3DISC.com