



Aidite®



Biomic™ Stain/Glaze

Simulation or mimic,
the final purpose is initial natural

■ Everyone with a healthy and beautiful smile

Biomic

High Aesthetics

Indications:

Zirconia Glass Ceramic

Sintering Chart:

Start temp.	Drying time	Heat rate	Start	End	Final temp.	Holding time	Slow cooling
450°C	8min	55°C/min	450°C	730°C	730°C	1min	-




Use Instruction

Accessories



 Simple to Use



Color Structure



Characteristic

High Saturation

Brown	Terracotta	Shade C	Shade D

Hue

Shade A	Shade B	Shade C	Shade D

Warm area

Orange 1	Yellow	Orange 2

Cold area

Mauve	Blue 1	Blue 2	Black


White	White-C	Brown




Main color area of gums




Intensive area of gums




Gum surface feature



Pink1, Pink2, Pink3



Pink1, Pink4



Pink5, 3D-T

 Bionic Aesthetics



Before Biomic



2D Effect




3D Effect




Final Effect

Fluorescent Glaze Paste



Effect of zirconia crown before
applying Glaze-F under fluorescent lamp



Effect of zirconia crown after
applying Glaze-F under fluorescent lamp

3D-Transpa & 3D-Value

3D paste can be used to adjust the tooth shape or contact area.

3D-Value can improve the brightness.

3D-Transpa can be mixed with 2D paste to create new 3D effect layer.



Glaze 3DT +



= No limit



A dental technician wearing glasses and a black shirt is working at a microscope in a laboratory setting. The technician is focused on a small object, likely a dental restoration, under the microscope's light. Various dental tools and materials are visible on the workbench.

Biomic

Stain/Glaze

Aesthetic
Gum

“

Biomic developed by Aidite completely allows me to finish the zirconia restorations without porcelain buildup, through the simplest and fastest way. Achieving a restoration just like the natural teeth.

Whether Biomic is applied to zirconia or glass ceramic restorations, it can achieve an incredibly beautiful and natural effect.

”

Kangni





 Biomic Stain/Glaze
Aesthetic



 Biomic Stain/Glaze
Gum

Aidite®

Aidite (Qinhuangdao) Technology Co.,Ltd.

Tel :0086-335-8587898

Web:www.aidite.com

Email:info@aidite.com

CFDA CE0197 FDAK111291



www.aidite.com